

BacLink

Data conversion and code mapping



**WHO Collaborating Centre for
Surveillance of Antimicrobial Resistance**

Boston, July 2022

Part 1. Running the conversion

You will now see how to convert your text file to a WHONET file with this new configuration. You will also see how to check whether the configuration is working and define any of your local data codes that WHONET doesn't recognize.

Starting the conversion


Original data file: Click on the File format for “BacLink Tutorial Hospital” that you configured in the previous step. In the upper box called “File name” click on “Browse” to select the file “ExcelDemo.txt”.



New data file


In the lower box, you need to indicate the name of the new WHONET file that you are about to create. For this tutorial, give the file the name: “ExcelDemo.bth”. You can give the new file any name that you like, but it will be convenient for your data file management if all of your data files end with the three-letter code that you gave to your laboratory, in this tutorial “bth”.

If you did these steps properly, your screen should look like the following.



Inspecting the conversion

Click on “Begin conversion”. BacLink will display for you results from the conversion of the first three isolates in the original data file. The purpose is to allow you to visually inspect the accuracy of the conversion. On the below screen, you see results from the first isolate. First focus on the middle column to see whether BacLink is reading the data values correctly, and check the final column to see whether BacLink is converting the data values correctly.




For this first isolate, there seems to be no problem with Identification number, Location, Date of birth, and Specimen number. There also appears to be no problem with Specimen date – you should always check the reformatting of the date to ensure that you have selected the correct date format (D/M/Y, M/D/Y, MMDDYYYY, Y-M-D, *etc.*). You will also notice that BacLink calculates the patient’s age using the date of birth and specimen date fields.


For Organism, you will notice that there are two similar rows: “Organism” and “Local organism code”. The first of these is intended for the WHONET organism code, while the second is intended to save the original organism code from your data file without modification. Since you have not yet defined the correspondence between your local code and the WHONET code, the first of these two rows is blank. The situation with Specimen type is similar. In the row “Specimen type”, BacLink leaves the WHONET field blank because it does not know what the corresponding WHONET code is for “Blood”. In the row “Local specimen code”, the original entry is copied over. Finally, in the lower box, you will see that BacLink was able to transfer the antibiotic results correctly.

So from this first record, we can conclude that BacLink is reading the data file correctly, and is saving most of the information into the new WHONET file accurately. However, BacLink has difficulty with the organism and specimen type fields because the codes have not yet been defined. If you think you may have made a mistake in your configuration, click “Cancel” and return to “Edit format” to make the correction. Otherwise, click on “Next” to proceed with the second isolate.

In the second isolate, you will notice a small problem that did not appear with the first isolate. The location code “Medicine 1” has been abridged to “medici”. This is because the default length of the location field is 6 characters. This will be easy to fix later by changing the field length to a larger value.




Click on “Next” to see the third isolate, and then “Next” again. BacLink will now proceed to finish the conversion of the rest of the data file. BacLink will tell you how many isolates it converted – 10 isolates in this tutorial. Click on “OK” to continue.




Defining unrecognized codes

When the conversion is finished, BacLink may alert you that it encountered some codes that it did not recognize. In this tutorial, you should see the below message. Click “Yes”.




BacLink will show you a summary of all of the different codes that it did not recognize, as in the below screen.




In this tutorial, BacLink indicates that there are unrecognized locations, organisms, and specimen types. We won't worry about the locations at this point until the problem with the field lengths is fixed (from six characters to a larger value).

Click on the row "Organism". Now click on "Define codes" to see the following screen with a list of all of the unrecognized organism codes/text.




Click on the first organism “C. albicans” and “Define code”. BacLink will now suggest a number of possible matches for this organism. If you see the correct match, click on the organism and click “OK”.




If you do not see the organism listed, then use the search box to look for the correct organism. For example, instead of “CNS” or “Coagulase-negative staphylococci”, WHONET has an entry “Staphylococcus, coagulase-negative”. To search for an organism, type one or more letters of the genus name and one or more letters of the species name, for example type “Staph coag” to find possible matches for “CNS”. For “GC”, you should do a search for something similar to “Neiss gon”.

After you match the first organism, continue to match all of the rest of the organisms. If you cannot find an appropriate match, you can call the organism “Other” or you can leave it undefined. Then click, “OK” to return to the list with the other undefined codes. You will see that the organism row has disappeared because all of the organisms have been defined.


Note: If you have a very long list of codes, there is no need to define all of them. Begin with the most common or important ones. If you do not define a code, WHONET will leave the corresponding WHONET field blank, but will still save the original organism code in the column “Local organism code”.



Click on the row for “Specimen type” and “Define codes”.



Define the specimen types one at a time, finding the closest match on the WHONET list. When you finish matching the specimen types, click “OK” to return to the remaining list of unrecognized codes.




Since we will define the locations later, click “Continue”. BacLink then warns that you have defined some additional codes that were not included during the conversion. Click “OK” to return to the main BacLink screen.



Running the conversion again...

Now that you have matched some of your codes with the appropriate WHONET codes, run the conversion again. Click on “Begin conversion”. BacLink will give you a warning (Warning #1) that the WHONET file that you want to create already exists, and asks whether you want to replace it. Answer “Yes”, since the first time was just an initial test.



BacLink will then ask again (Warning #2) just to make sure that you read the question. Answer “Yes” again.

You will now see the first isolate again, but this time with the appropriate WHONET organism and specimen type codes indicated. So this conversion is more complete than during your first attempt.

Field name	ExcelDemo.txt	ExcelDemo.bth
► Identification number	1658752	1658752
Last name		
First name		
Full name		
Sex		
Date of birth	3/10/1957	03-Oct-1957
Age		48
Location	Ortho	ortho
Department		
Specimen number	45896	45896
Specimen date	22/12/2005	22-Dec-2005
Specimen type	Blood	bl
Local specimen code	Blood	Blood
Organism	K_pneumo	kpn
Local organism code	K_pneumo	K_pneumo

Antibiotics


AMP_ND10	6
CIP_ND5	23
SXT_ND1_2	24

Next Cancel


Click “Next isolate”, and you will see that the problem with the length of the location code field still persists (“medici” instead of “medicine 1”). Click “Next” and “Next” again to finish the conversion. BacLink indicates that it converted ten isolates. Click “OK”, and answer “Yes” to the question about the unrecognized codes. This time, you will see that the list of unrecognized codes is much shorter than before. Only the locations remain to be defined. Click on “Continue” to return to the main BacLink screen.

Editing the configuration

The conversion is working in most respects with the exception of the location codes which are being shortened to six characters. This is easy to fix. Click on “Edit format” from the main BacLink screen to return to the configuration part of BacLink. Click on “File format” and “Data fields”. This will return you to the screen in which you defined the matches between the fields in the Excel data file and in WHONET.




Click on “Modify the list of data fields”. You will see the following screen. From this screen, you can configure the current fields, add additional ones to the list, and delete fields that you do not need. Click on the option “Location”. On the right side of the screen, you will see that this length of this field is set to 6 characters. Change this to 20.




Then click “OK”, “OK”, “OK”, “Save”, and “Exit”. This should return you to the main BacLink screen.

Running the conversion again...

Click on “Begin conversion”. Answer “Yes” and “Yes” again to replace the current WHONET file. BacLink will show you the first isolate, but this is unchanged from before. Click on “Next” to see the second isolate. You will now notice that the location “Medicine 1” is being correctly being saved as “medicine 1”.



Then continue with the rest of the conversion until you see the list of unrecognized codes. Click on “Location” and “Define codes”.



For the first item on the list, click on “Define code” to get the following screen.

You will notice that defining location codes is different from defining organism, specimen, and antibiotic codes. For the latter three, WHONET has a standard list of codes used by all laboratories with WHONET. For location codes, WHONET simply uses the same code that your laboratory uses.

Defining your location codes serves a different purposes – it allows you to describe the details about the location to BacLink, for example which department/specialty, whether the location is in an inpatient or outpatient area, *etc.*

For the first location, “5 North”, choose a department of “med=Medicine” and a location type of “inx = Inpatient (non-ICU)”. Then click “OK”.

Continue defining the other location codes in the manner using the following suggestions:

<u>Location</u>	<u>Institution</u>	<u>Department</u>	<u>Location type</u>
Diabetes clinic	bth	med	out
ER	bth	eme	eme
HC5	oth	out	out
ICU	bth	icu	icu
Medicine 1	bth	med	inx
Medicine 2	bth	med	inx
NICU	bth	neo	icu
OB-Gyn	bth	obg	inx

Ortho

bth

sur

in


In this example, HC5 represents a “Health Center #5”, an outpatient clinic outside the hospital.

Defining locations has two important benefits: 1. at the local level, hospital staff will be able to analyze easily data at the level of the ward, department, or inpatient/outpatient status, even if the original data file only had the ward information; and 2. at the national level, defining the location codes permits standardization for purposes of national comparisons and analyses.

After defining all of the location codes, return to the main BacLink screen.

Running the conversion again... for the last time!

Click on “Begin conversion”, and replace the previous WHONET file. Look at the conversion for the first isolate. You will notice that for the location “Ortho”, BacLink is now putting the department equal to “sur=Surgey” because of the location matchings you defined in the previous step.



Field name	ExcelDemo.txt	ExcelDemo.bth
Identification number	1658752	1658752
Last name		
First name		
Full name		
Sex		
Date of birth	3/10/1957	03-Oct-1957
Age		48
Location	Ortho	ortho
Department		sur
Specimen number	45896	45896
Specimen date	22/12/2005	22-Dec-2005
Specimen type	Blood	bl
Local specimen code	Blood	Blood
Organism	K. pneumo	kpn
Local organism code	K. pneumo	K. pneumo

Antibiotics

AMP_ND10	6
CIP_ND5	23
SXT_ND1_2	24


Then click “Next” and finish the rest of the conversion. When the conversion is finished, BacLink will return you immediately to the main BacLink screen because all of the undefined codes have now been defined.

Congratulations! You have successfully finished a BacLink configuration, defined all of the codes needed by WHONET, and converted a sample data file (ExcelDemo.xls) to WHONET format (ExcelDemo.bth). You can then proceed to WHONET to analyze this file. Click on “Exit” to leave BacLink.

The configuration that you made will be useful for any file that has this same structure. For example, if use your January data file to define your BacLink configuration and codes matches, you can use the same configuration to convert your February, March, *etc.* data. If BacLink encounters some new codes in the following months, you can define them as they arise over time.

Part 2. Getting started with WHONET

Now that you have created a valid WHONET file using BacLink, you can proceed to WHONET. Double-click on your WHONET icon to begin WHONET.




You will notice that the “BacLink Tutorial Hospital” that we are using in this tutorial does not yet appear on the WHONET list of laboratories. Before you can begin analyzing the sample data file, you will need to create a laboratory configuration


Creating a laboratory configuration

To create a WHONET laboratory configuration, one option would be to use the “New laboratory” option on this screen. You would need to enter the laboratory name and code, the list of antibiotics, locations, *etc.* But fortunately, there is a faster way.


Click on “Cancel” to go to the main WHONET menu. Click on “File”. You will notice an option called “Create a laboratory from a data file”. Click on this option.




Enter the country “World Health Organization”, laboratory = “BacLink Tutorial Hospital”, and laboratory code = “BTH”. Click “OK”.




On the next screen, WHONET will ask you to select an example of a data file from this laboratory. Choose the file “ExcelDemo.bth” that you created with this BacLink tutorial, and click “Open”.




WHONET will now examine the file “ExcelDemo.bth” – what antibiotics are present, what data fields, locations, *etc.* With this information WHONET can create a matching laboratory configuration. When WHONET finishes, the program will ask if you want to review the new configuration. Answer “Yes”.




WHONET will then show you the laboratory configuration that it prepared.




If you click on “Antibiotics”, you will see that WHONET has the three antibiotics present in the original Excel file. You can also check the antibiotic breakpoints, and configure the panels and profiles if you wish.



If you click on “Locations”, you will see the list of locations from the original file.



If you click on “Data fields”, you will see that the location field has a length of 20 characters, rather than the usual default of 6 characters.




After you explore these screens, click on “Save” to leave the configuration program.

Note: After creating the configuration utilizing the shortcut described here, you can make any further edits, such as any modifications to the antibiotic breakpoints and profiles, with **Modify laboratory** in the same way as any WHONET laboratory configuration.

Inspecting the data file

As a final step, click on “Data entry”, “Open data file”, and choose the file “ExcelDemo.bth”. You will now be in the Data entry program. Click on “View database” to see the following screen.




The screenshot shows a window titled "Data entry: C:\whonet5\Data\ExcelDemo.bth". The window contains a table with the following columns: Identification number, Specimen number, Country, Laboratory, Origin, Last name, First name, Sex, Date of birth, Age, Location, Institution, and Department. The data is as follows:

Identification number	Specimen number	Country	Laboratory	Origin	Last name	First name	Sex	Date of birth	Age	Location	Institution	Dept
1698752	45896	WHO	BTH					10/3/1957	48	ortho	bth	sur
4693314	32047	WHO	BTH					3/15/1980	25	medicine1	bth	med
7584301	14658	WHO	BTH					2/2/1940	65	medicine2	bth	med
4865921	76481	WHO	BTH					1/10/2006		diabetes clinic	bth	med
5481023	16594	WHO	BTH					6/20/1924	81	icu	bth	icu
4935814	68459	WHO	BTH					2/6/1992	13	ob-gyn	bth	obg
9208748	28791	WHO	BTH					7/26/2004		5 north	bth	med
3017945	30247	WHO	BTH					11/30/2001	4	hc5	oth	out
4869704	45076	WHO	BTH					9/11/1952	53	er	bth	eme
5904720	48652	WHO	BTH					12/18/1966	19	nicu	bth	neo

WHONET is now showing you the complete contents of the data file created by BacLink. Inspect this file to see if there are any possible errors – in the dates, codes, antibiotic results, *etc.* This will also help to familiarize you with the structure and contents of a WHONET data file.

When you finish, click on “Continue” and “Exit” to return to the main WHONET screen. Then “File”, “Exit” to leave WHONET completely.

The next time that you enter WHONET, you will see your new laboratory configuration for the “BacLink Tutorial Hospital” on the list. You can then select this laboratory for laboratory configuration, data entry, or data analysis.



Part 3. What the next steps?

You have learned how to convert data from common desktop softwares, such as Excel, and text files exported from laboratory information systems. The steps were: 1. configuration – which only needs to be accomplished once; and 2. running the conversion – which you can repeat as you continue to collect additional data.

You then set up a new WHONET configuration for this laboratory. With this configuration, you can now analyze your data with WHONET. If needed, you could also use WHONET for manual data entry of supplemental information that you did not have in your original data file.

For guidance on how to use WHONET for data entry and data analysis, continue with the appropriate tutorials.